

# SIC presents. Optical Waveguide spectroscopy

## SIS5000/S-SPR6000 series

**Optical waveguide spectrometry** is very sensitive, since the incident light advances by repeated total reflection in the waveguide substrate, and a sample on the waveguide surface absorbs the **evanescent wave** at every reflection point. The intensity of the **evanescent wave** diminishes strongly with distance from waveguide surface, and only distance similar to the wavelength of the incident light are effective in detecting the spectrum. This **high-sensitivity** is one of advantage for **interfacial analysis**.




Fig.2 sample room for analysis




Fig.1 Image of optical waveguide spectrometry




Fig.4 wave-guide variation




Fig.3 S-SPR6000 system front view

### Applications

- Photo functional surface, Photo chromism, Photo deplete
- Collaborate electrochemistry, collaborate fluorescence spectra
- Molecule devices, Organic EL, Dye-sensitized Solar Cell
- Bio sensing by spectroscopic-SPR, interaction of bio molecule
- Langmuir-Blodgett membrane, molecule orientation ...etc



## ◇ Specification

### ■ The main system

Incident angle adjustment	normal to 90~35degrees(fine adjustment bellow 0.005)
Out put angle adjustment	normal to 90~35degrees(fine adjustment bellow 0.005)
Indication choose	normal to angle or parallel to angle
Y scan adjustment	0~20mm(fine adjustment bellow 0.05)
Light source	Xe 150W(Hamamatsu photonics KK L2175)
Dimension	620(W)×550(D)×330(H)

### ■ Computer

Windows XP Base (Dell computer Dimension well)  
Original application software for absorption spectra and SPR spectra and save on time

■ Power 100~110V AC 0.6KVA

■ Spectrometer Ocean Optics inc S2000(S4000) spectrometer

Wavelength 220~700nm or 400~1000nm ,  
220~700nm and 400~1000nm 2 spectrometers


CCD 2048pix , resolution 1.25nm

### ■ Optical wave-guide type

Dimension 0.2(t)×65(l)×20(W) each side has 60 degree tilt angle.

Quartz or high reflective index glass attached Au, attached ITO

### ■ Options



Transmission cell




Z axis adjust stage




Polarizer unit



Halogen light source



electro chemical cell



PDMS μ-line

Manufacturer

**SIC** SYSTEM INSTRUMENTS CO., LTD.



**EVER SEIKO CORPORATION**

ES Bldg., 4-39-5, Senzoku, Taitoku, Tokyo 111-0031 Japan

Tel: (813) 3874-8871

URL: <http://www.everseiko.co.jp>

Fax: (813) 3874-8873

E-mail: [info@everseiko.co.jp](mailto:info@everseiko.co.jp)